

**AZENTA**  
LIFE SCIENCES

azenta.com


# BIOINFORMATICS REPORT

Pacbio De-novo Hifi Eukaryotic Assembly

**AZENTA LIFE SCIENCES**  
**NEXT GENERATION SEQUENCING**


115 Corporate Boulevard  
South Plainfield, NJ 07080  
Email: [ngs@azenta.com](mailto:ngs@azenta.com)  
Phone: 908-222-0711 Ext. 1

5/23/2024



## Table of Contents

|          |   |          |
|----------|---|----------|
| <b>1</b> | <b>Project Information</b>                  | <b>2</b> |
| <b>2</b> | <b>Description of Workflow</b>              | <b>3</b> |
| 2.1      | Description of library preparation workflow | 3        |
| 2.2      | Workflow of Data Analysis                   | 4        |
| <b>3</b> | <b>Data Analysis Results</b>                | <b>5</b> |
| 3.1      | Assembly Statistics                         | 6        |
| <b>4</b> | <b>Deliverables</b>                         | <b>7</b> |



**AZENTA**  
LIFE SCIENCES

## 1 Project Information


|             |  |
|-------------|--|
| Project     | 30-742644819   |
| Customer    | Thomas Parchman  |
| Institution | University of Nevada - Reno                              |
| Email       | <a href="mailto:tparchman@unr.edu">tparchman@unr.edu</a> |
| Service     | Pacbio Hifi DNA De-novo Assembly Package                 |
| Species     | <i>Loxia curvirostra</i>                                 |
| Sample      | 2  |




## 2 Description of Workflow

### 2.1 Description of library preparation workflow


PacBio library was prepared using SMRTbell prep kit per the manufacturer's protocol. Follow the run design guidelines for sequencing. Figure 1 outlines the PacBio library preparation and sequencing workflow.



*Figure 1: Pacbio library preparation and sequencing workflow 3. CCS ANALYSIS*




## 2.2 Workflow of Data Analysis





### 3 Data Analysis Results

The Hifi data for each of the samples was de-novo assembled to a genome using FlyE v2.9.3-b1797. The raw data was assessed for quality trimming using fastp v0.23.4 setting of length being 2000 bps and quality score of each read being a minimum of 15. FlyE v2.9.3-b1797 was used with default settings with polishing to generate a draft de-novo assembly. We then filtered the draft assembly for small contigs (less than 50Kbps) to generate a final assembly. The final assembly was analyzed for open reading frames via EMBOSS tools ORF Finder. The identified protein ORF was used for Diamond BLASTp annotation through NCBI's NR database. Also, Quast was then used to generate general assembly statistics.



### 3.1 Assembly Statistics

Quast Statistics for each assembly.

| <b>Assembly Statistics</b>      | 171203-5      | S170996-2     |
|---------------------------------|---------------|---------------|
| # contigs ( $\geq 0$ bp)        | 9,765         | 9,014         |
| # contigs ( $\geq 1000$ bp)     | 9,765         | 9,014         |
| # contigs ( $\geq 5000$ bp)     | 9,765         | 9,014         |
| # contigs ( $\geq 10000$ bp)    | 9,765         | 9,014         |
| # contigs ( $\geq 25000$ bp)    | 9,765         | 9,014         |
| # contigs ( $\geq 50000$ bp)    | 9,765         | 9,014         |
| Total length ( $\geq 0$ bp)     | 1,495,363,063 | 1,697,079,343 |
| Total length ( $\geq 1000$ bp)  | 1,495,363,063 | 1,697,079,343 |
| Total length ( $\geq 5000$ bp)  | 1,495,363,063 | 1,697,079,343 |
| Total length ( $\geq 10000$ bp) | 1,495,363,063 | 1,697,079,343 |
| Total length ( $\geq 25000$ bp) | 1,495,363,063 | 1,697,079,343 |
| Total length ( $\geq 50000$ bp) | 1,495,363,063 | 1,697,079,343 |
| # contigs                       | 9,765         | 9,014         |
| Largest contig                  | 3,863,186     | 5,613,873     |
| Total length                    | 1,495,363,063 | 1,697,079,343 |
| GC (%)                          | 42.64         | 43.38         |
| N50                             | 197,710       | 286,440       |
| N75                             | 101,361       | 130,809       |
| L50                             | 1,817         | 1,439         |
| L75                             | 4,493         | 3,698         |
| # N's per 100 kbp               | 0.85          | 0.72          |



## 4 Deliverables

- Standard HiFi Sequel IIe
  - Hifireads.bam
  - Hifireads.fastq.gz
- Assembly
  - FlyE v2.9.3-b1797 – filtered for small contigs (50 Kbps)
    - Assembled contigs (fasta)
  - FlyE v2.9.3-b1797 - Unfiltered includes small contigs (50 Kbps)
    - Unfiltered contigs (fasta)
- Annotation
  - EMBOSS ORF Finder
    - Nucleotide ORFs – fasta and table (TSV)
    - Protein ORFs – fasta and table (TSV)
  - Diamond BLASTp
    - Annotation Table (TSV)
- Statistics
  - FlyE v2.9.3-b1797
    - Per contig stats on unfiltered contigs (TXT)
  - Quast – Assembly statistics
    - General assembly statistics on filtered de-novo assembly (TSV)